

## Product Datasheet

### Human C-MAPT/C-TAU ELISA Kit orb404157

<b>Description</b>	Human cleaved microtubule-associated protein tau (C-MAPT/C-TAU) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	62.5 pg/mL-4000 pg/mL
<b>Sensitivity</b>	15.6 pg/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	C-MAPT/C-TAU
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, cerebrospinal fluid (CSF), tissue homogenates, cell lysates
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	• Sample Type: serum, plasma, cerebrospinal fluid (CSF), tissue homogenates, cell lysates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450nm