

Product Datasheet

Human ST3GAL1 ELISA Kit [Out of Stock] orb404159

Description	Human CMP-N-acetylneuraminate-beta-galactosamide-alpha-2,3-sialyltransferase 1(ST3GAL1) ELISA kit
Reactivity	Human
Range	18.75 pg/ml-1200 pg/ml
Sensitivity	4.67 pg/ml.
Tested Applications	ELISA
Target	ST3GAL1
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates, cell lysates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm