

Product Datasheet

Human COL8A2 ELISA Kit orb404186

Description	Human Collagen alpha-2(VIII) chain (COL8A2) ELISA kit
Reactivity	Human
Range	0.625 ng/ml-40 ng/ml
Sensitivity	0.156 ng/ml
Tested Applications	ELISA
Target	COL8A2
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm