

Product Datasheet

Human Collagen I ELISA Kit orb404196

Description	Human Collagen I ELISA Kit
Reactivity	Human
Range	0.625 ng/ml-40 ng/ml
Sensitivity	
Tested Applications	ELISA
Target	Collagen I
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, cell culture supernates, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm