

Product Datasheet

Human Complement component C6 ELISA Kit orb404210

Description	Human Complement component C6(C6) ELISA kit
Reactivity	Human
Range	62.5 pg/mL-4000 pg/mL
Sensitivity	15.6 pg/mL
Tested Applications	ELISA
Target	C6
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P13671
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm