

Product Datasheet

Human CLEC9A ELISA Kit orb404236

Description	Human C-type lectin domain family 9 member A(CLEC9A) ELISA kit
Reactivity	Human
Range	0.312 ng/ml-20 ng/ml
Sensitivity	0.078 ng/ml.
Tested Applications	ELISA
Target	CLEC9A
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates, cell lysates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm