

Product Datasheet

Human DPEP1 ELISA Kit orb404300

Description	Human Dipeptidase 1(DPEP1) ELISA kit
Reactivity	Human
Range	23.5 pg/mL-1500 pg/mL
Sensitivity	5.86 pg/mL
Tested Applications	ELISA
Target	DPEP1
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P16444
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates, cell lysates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm