

## Product Datasheet

### Human CD248 ELISA Kit orb404367

<b>Description</b>	Human Endosialin(CD248) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	0.312 ng/mL-20 ng/mL
<b>Sensitivity</b>	0.078 ng/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	CD248
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates
<b>Uniprot ID</b>	<b>Q9HCU0</b>
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	• Sample Type: serum, plasma, tissue homogenates, cell lysates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm