

## Product Datasheet

### Human EGFL6 ELISA Kit orb404380

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Human Epidermal growth factor-like protein 6(EGFL6) ELISA kit   |
| <b>Reactivity</b>          | Human   |
| <b>Range</b>               | 78 pg/ml-5000 pg/ml   |
| <b>Sensitivity</b>         | 19.5 pg/ml.   |
| <b>Tested Applications</b> | ELISA   |
| <b>Target</b>              | EGFL6   |
| <b>Note</b>                | For research use only.  |
| <b>Expiration Date</b>     | Please enquire.   |
| <b>Application Notes</b>   | <ul style="list-style-type: none"><li>• <b>Sample Type:</b> serum, plasma, tissue homogenates</li><li>• <b>Assay Time:</b> 1-5 hours</li><li>• <b>Sample Volume:</b> 50-100 <math>\mu</math>l</li><li>• <b>Detection Wavelength:</b> 450 nm</li></ul> |