

Product Datasheet

Human FAM12B ELISA Kit orb404384

Description	Human Epididymal secretory protein E3-beta(FAM12B) ELISA kit
Reactivity	Human
Range	0.1 ng/ml-20 ng/ml
Sensitivity	0.05 ng/ml
Tested Applications	ELISA
Target	FAM12B
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm