
Product Datasheet

Human EPGN ELISA Kit (orb404390)

Description	Human Epigen(EPGN) ELISA kit
Reactivity	Human
Range	1.56 ng/mL-100 ng/mL
Sensitivity	0.39 ng/mL
Tested Applications	ELISA
Target	EPGN
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q6UW88
Expiration Date	Please enquire.
Application Notes	• Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm