

Product Datasheet

Human Glucocorticoid Receptor ELISA Kit orb404484

Description	Human Glucocorticoid receptor(NR3C1) ELISA kit
Reactivity	Human
Range	31.25 pg/ml-2000 pg/ml
Sensitivity	7.81 pg/ml
Tested Applications	ELISA
Target	Glucocorticoid Receptor
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates, cell lysates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm