

Product Datasheet

Human GH2 ELISA Kit orb404521

Description	Human Growth hormone variant(GH2) ELISA kit
Reactivity	Human
Range	133 pg/ml -4000 pg/ml
Sensitivity	33pg/ml
Tested Applications	ELISA
Target	GH2
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm