

## Product Datasheet

### Human HDAC9 ELISA Kit orb404564

<b>Description</b>	Human Histone deacetylase 9(HDAC9) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	0.156 ng/ml-10 ng/ml
<b>Sensitivity</b>	0.039 ng/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	HDAC9
<b>Note</b>	For research use only.
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	<ul style="list-style-type: none"><li>• <b>Sample Type:</b> serum, plasma, tissue homogenates, cell lysates</li><li>• <b>Assay Time:</b> 1-5 hours</li><li>• <b>Sample Volume:</b> 50-100 <math>\mu</math>l</li><li>• <b>Detection Wavelength:</b> 450 nm</li></ul>