

Product Datasheet

Human CD74 ELISA Kit orb404573

Description	Human HLA class II histocompatibility antigen gamma chain(CD74) ELISA kit
Reactivity	Human
Range	0.9 ng/mL-60 ng/mL
Sensitivity	0.225 ng/mL
Tested Applications	ELISA
Target	CD74
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P04233
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm