

## Product Datasheet

### Human INSL5 ELISA Kit orb404618

<b>Description</b>	Human Insulin-like peptide INSL5(INSL5) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	8 ng/mL-140 ng/mL
<b>Sensitivity</b>	5 ng/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	INSL5
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>Q9Y5Q6</b>
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	• Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm