

Product Datasheet

Human ITGB8 ELISA Kit (orb404627)

Description	Human Integrin beta-8(ITGB8) ELISA kit
Reactivity	Human
Range	0.39 ng/mL-25 ng/mL
Sensitivity	0.09 ng/mL
Tested Applications	ELISA
Target	ITGB8
Note	For research use only.
Sample Types	serum, plasma,tissue homogenates
Uniprot ID	P26012
Expiration Date	Please enquire.
Application Notes	• Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm