

Product Datasheet

Human IL17RB ELISA Kit orb404655

Description	Human Interleukin-17 receptor B(IL17RB) ELISA kit
Reactivity	Human
Tested Applications	ELISA
Target	IL17RB
Note	For research use only.
Uniprot ID	Q9NRM6
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates, cell lysates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm