

Product Datasheet

Human IL3RA ELISA Kit orb404669

Description	Human Interleukin-3 receptor subunit alpha(IL3RA) ELISA kit
Reactivity	Human
Range	0.312 ng/ml-20 ng/ml
Sensitivity	0.078 ng/ml.
Tested Applications	ELISA
Target	IL3RA
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates, cell lysates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm