

Product Datasheet

Human IL1F5 ELISA Kit orb404673

Description	Human Interleukin-36 receptor antagonist protein(IL1F5) ELISA kit
Reactivity	Human
Tested Applications	ELISA
Target	IL1F5
Note	For research use only.
Uniprot ID	Q9UBH0
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, cell culture supernates, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm