

Product Datasheet

Human KEAP1 ELISA Kit orb404686

Description	Human Kelch-like ECH-associated protein 1(KEAP1) ELISA kit
Reactivity	Human
Range	25 pg/ml-1600 pg/ml
Sensitivity	6.25 pg/ml
Tested Applications	ELISA
Target	KEAP1
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates, cell lysates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm