
Product Datasheet

Human LAMC1 ELISA Kit (orb404705)

Description	Human Laminin subunit gamma-1(LAMC1) ELISA kit
Reactivity	Human
Range	0.78 ng/mL-50 ng/mL
Sensitivity	0.195 ng/mL
Tested Applications	ELISA
Target	LAMC1
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P11047
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm