

Product Datasheet

Human LGR5 ELISA Kit orb404724

Description	Human Leucine-rich repeat-containing G-protein coupled receptor 5(LGR5) ELISA kit
Reactivity	Human
Sensitivity	4.69 pg/ml.
Tested Applications	ELISA
Target	LGR5
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm