

## Product Datasheet

### Human LILRA3 ELISA Kit (orb404733)

<b>Description</b>	Human Leukocyte immunoglobulin-like receptor subfamily A member 3(LILRA3) ELISA kit...
<b>Reactivity</b>	Human
<b>Range</b>	0.156 ng/mL-10 ng/mL
<b>Sensitivity</b>	0.039 ng/mL
<b>Target</b>	LILRA3
<b>Note</b>	For research use only.
<b>Sample Types</b>	serum, plasma, tissue homogenates
<b>Uniprot ID</b>	<b>Q8N6C8</b>
<b>Expiration Date</b>	Please enquire.