

Product Datasheet

Human FCGR3A ELISA Kit (orb404761)

Description	Human Low affinity immunoglobulin gamma Fc region receptor III-A(FCGR3A) ELISA kit
Reactivity	Human
Range	47 pg/mL-3000 pg/mL
Sensitivity	11.75 pg/mL
Tested Applications	ELISA
Target	FCGR3A
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P08637
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm