
Product Datasheet

Human LYPD1 ELISA Kit (orb404767)

Description	Human Ly6/PLAUR domain-containing protein 1(LYPD1) ELISA kit
Reactivity	Human
Range	25 pg/mL-1600 pg/mL
Sensitivity	6.25 pg/mL
Tested Applications	ELISA
Target	LYPD1
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q8N2G4
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates, cell lysates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm