

Product Datasheet

Human LIPA ELISA Kit (orb404777)

Description	Human Lysosomal acid lipase/cholesteryl ester hydrolase(LIPA) ELISA kit
Reactivity	Human
Range	62.5 pg/mL-4000 pg/mL
Sensitivity	15.6 pg/mL
Tested Applications	ELISA
Target	LIPA
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P38571
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates and cell lysates• Assay Time: 1-5 hours• Sample Volume: 50-100 µl• Detection Wavelength: 450 nm