

Product Datasheet

Human SCGB2A1 ELISA Kit orb404793

Description	Human Mammaglobin-B(SCGB2A1) ELISA kit
Reactivity	Human
Range	0.235 ng/mL-15 ng/mL
Sensitivity	0.058 ng/mL
Tested Applications	ELISA
Target	SCGB2A1
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O75556
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm