

Product Datasheet

Human MUSK ELISA Kit orb404874

Description	Human Muscle Skeletal Receptor Tyrosine Kinase (MUSK) ELISA Kit
Reactivity	Human
Range	1.56 pg/ml -100 pg/ml
Sensitivity	0.39pg/ml
Tested Applications	ELISA
Target	MUSK
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm