

Product Datasheet

Human NLRP1 ELISA Kit orb404892

Description	Human NACHT, LRR and PYD domains-containing protein 1(NLRP1) ELISA kit
Reactivity	Human
Range	18.75 pg/ml-1200 pg/ml
Sensitivity	4.67 pg/ml
Tested Applications	ELISA
Target	NLRP1
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates, cell lysates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm