

Product Datasheet

Human MT-ND6 ELISA Kit orb404895

Description	Human NADH-ubiquinone oxidoreductase chain 6(MT-ND6) ELISA kit
Reactivity	Human
Range	23.44 pg/ml-1500 pg/ml
Sensitivity	5.86 pg/ml
Tested Applications	ELISA
Target	MT-ND6
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates, cell lysates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm