

Product Datasheet

Human NAPSA ELISA Kit orb404898

Description	Human Napsin-A(NAPSA) ELISA kit
Reactivity	Human
Range	4.7 ng/mL-300 ng/mL
Sensitivity	1.18 ng/mL
Tested Applications	ELISA
Target	NAPSA
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O96009
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm