

Product Datasheet

Human DCC ELISA Kit [Out of Stock] orb404910

Description	Human Netrin receptor DCC(DCC) ELISA kit
Reactivity	Human
Range	47 pg/ml-3000 pg/ml
Sensitivity	11.75 pg/ml
Tested Applications	ELISA
Target	DCC
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates and cell lysates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm