

Product Datasheet

Human RAET1E ELISA Kit [Out of Stock] orb404940

Description	Human NKG2D ligand 4(RAET1E) ELISA kit
Reactivity	Human
Range	1.56 ng/ml-100 ng/ml
Sensitivity	0.39ng/ml
Tested Applications	ELISA
Target	RAET1E
Note	For research use only.
Expiration Date	Please enquire.
Application Notes	<ul style="list-style-type: none">• Sample Type: serum, plasma, tissue homogenates• Assay Time: 1-5 hours• Sample Volume: 50-100 μl• Detection Wavelength: 450 nm