

## Product Datasheet

### Human NOD1 ELISA Kit orb404964

<b>Description</b>	Human Nucleotide-binding oligomerization domain-containing protein 1(NOD1) ELISA kit
<b>Reactivity</b>	Human
<b>Range</b>	15.6 pg/ml-1000 pg/ml
<b>Sensitivity</b>	3.9pg/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	NOD1
<b>Note</b>	For research use only.
<b>Expiration Date</b>	Please enquire.
<b>Application Notes</b>	<ul style="list-style-type: none"><li>• <b>Sample Type:</b> serum, plasma, tissue homogenates, cell lysates</li><li>• <b>Assay Time:</b> 1-5 hours</li><li>• <b>Sample Volume:</b> 50-100 <math>\mu</math>l</li><li>• <b>Detection Wavelength:</b> 450 nm</li></ul>