

Product Datasheet

Human PPIC ELISA Kit (orb404991)

Description	Human Peptidyl-prolyl cis-trans isomerase C(PPIC) ELISA kit
Reactivity	Human
Range	23.5 pg/mL-1500 pg/mL
Sensitivity	5.8 pg/mL
Tested Applications	ELISA
Target	PPIC
Note	For research use only.
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P45877
Expiration Date	Please enquire.
Application Notes	• Sample Type: serum, plasma, tissue homogenates • Assay Time: 1-5 hours • Sample Volume: 50-100 µl • Detection Wavelength: 450 nm